


Delfi Pro back parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	285299
2	158810
3	144433
4	111014
5	285261 (right, size 1) 285263 (right, size 2) 285265 (right, size 3)
6	167625
7	285254 (size 1-2) 285255 (size 3)
8	285262 (left, size 1) 285264 (left, size 2) 285266 (left, size 3)
9	124416
10	285257 (size 1) 285258 (size 2) 285259 (size 3)
11	285250 (size 1) 285251 (size 2) 285252 (size 3)
12	152665

Delfi Pro leg parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	127626
2	129620
3	129620
4	135009
5	285241 (left, size 1) 285243 (left, size 2) 285245 (left, size 3) 285240 (right, size 1) 285242 (right, size 2) 285244 (right, size 3)
6	285207 (left, size 1) 285210 (left, size 2) 285213 (left, size 3) 285206 (right, size 1) 285209 (right, size 2) 285212 (right, size 3)

Delfi Pro seating parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
7	285246 (size 1) 285247 (size 2) 285248 (Size 3)
8	285375
9	285294 (size 1) 285295 (size 2) 285296 (size 3)
10	285375
11	285308 (size 1) 285311 (size 2) 285314 (size 3)
12	285233
13	125204
14	285351 (size 1) 285352 (size 2) 285353 (size 3)

Delfi Pro seating base


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	166625
2	285225
3	135008
4	285227
5	915000 (size 1) 915022 (size 2)
6	285231
7	285229
8	135006
9	130812
10	285232
11	285236 (size 1) 285237 (size 2) 285238 (size 3)
12	127612
13	285221 (size 1) 285222 (size 2) 285223 (size 3)
14	285234
15	166614
16	135008
17	130812
18	285226
19	285230

Delfi Pro headrest parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	283016 (size 1) 283017 (size 1, long) 283018 (size 2) 283019 (size 2, long)
2	285381 (size 1) 285382 (size 1, long) 285383 (size 2) 285384 (size 2, long)
3	283110
4	152210
5	159515

Delfi Pro armrest parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	129612
2	142630
3	141612
4	990065
5	285276 (size 1) 285277 (size 2) 285278 (size 3)
6	162005
7	160600
8	158810
9	285275 (size 1) 285279 (size 2-3)


WHEN ORDERING THE ARMREST PAD FOR OLDER ARMREST BRACKETS

The measure A (c/c) on your existing bracket should be 50 mm.
If not you need to also order this bracket.



Delfi Pro footrest parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro



Item	Article number
1	27525 (size 1) 27526 (size 2) 27527 (size 3)
2	127626
4	166614
5	144433
6	285309 (size 1) 285310 (size 2-3)
7	285320
8	111005
9	152613
10	285317 (right) 285318 (left)

Delfi Pro footrest parts

Anatomic SITT AB Exploded views of Zitzi Delfi Pro




Item	Article number
1	213130
2	213110
3	285371 (size 1) 285372 (size 2) 285373 (size 3)

Delfi Pro footrest parts

Anatomic SITT AB Exploded views of Zitzi Delfi Pro


Item	Article number
1	152291
2	213170
3	159510



Delfi Pro footrest parts

Anatomic SITT AB Exploded views of Zitzi Delfi Pro


Item	Article number
1	285361
2	285361
3	213150



Delfi Pro footrest parts


Anatomic SITT AB Exploded views of Zitzi Delfi Pro

Item	Article number
1	152456



Delfi Pro calf support

Anatomic SITT AB Exploded views of Zitzli Delfi Pro



Item	Article number
1	283108 (size 1-2) 283109 (size 3)
2	283111
3	152613
4	127612
5	166614
6	135008

Overview of mounting- & adapter systems

Adapter system	fits on frame	with seating system	The seating system is
285171	Seabass Pro / Flipper Pro	Delfi Pro size 1	Removable
285172	Seabass Pro / Flipper Pro	Delfi Pro size 2	Removable
285173	Seabass Pro / Flipper Pro	Delfi Pro size 3	Removable
285174	R82	Delfi Pro size 1-3	Removable
285175 ¹⁾	Universal system	Delfi Pro size 1-3	Removable
285199	R82	Delfi Pro size 1-3	Removable
285199	R82	Sharky Pro size 1-2	Removable
153040 ²⁾	Universal system	Sharky Pro size 1-2	Removable
75133	Early Anatomic SITT	Delfi size 1, Anatomic Sits	Removable
75134	Early Anatomic SITT	Delfi size 2, Anatomic Sits	Removable
75135	Early Anatomic SITT	Delfi size 3, Anatomic Sits	Removable
26531 ³⁾	Universal system	Delfi size 1, Anatomic Sits	Removable
26532 ³⁾	Universal system	Delfi size 2, Anatomic Sits	Removable
26533 ³⁾	Universal system	Delfi size 3, Anatomic Sits	Removable
75137	R82	Delfi size 1, Anatomic Sits	Removable
75138	R82	Delfi size 2, Anatomic Sits	Removable
75139	R82	Delfi size 3, Anatomic Sits	Removable

















1) Contains of 152186 & 285176

2) Contains of 152186 & 153041

3) Strengthened with a wooden bottom plate

Removable means that you can separate the seating unit and the frame from each other without using of any tools.





Tip!

If you want to move the seating system between multiple chassis - buy an extra 152186 when usage of 285176/153041.

