

"Flex-joint"

The seat is equipped with a "Flex-joint" to make the sitting position more dynamic, and this helps to break any active extensor tone and/or spasticity. The angle of the back can be adjusted from 50 forwards to 30 backwards (starting position 90° hip angle).



Adjustable back

This is recommended for clients with a scoliosis. The thoracic pieces are independently adjustable in height and width.



Adjustable leg pieces

The leg pieces can be adjusted to accommodate independent hip abduction or adduction, as well as specific leg length discrepancies. The leg pieces can also be adjusted to accommodate a fixed pelvic obliquity.



Infinite abduction



Leg length difference