





## TECHNICAL DATA

Zitzi Delfi	Size 1	Size 2	Size 3
User weight, approx.:	25 kg	30 kg	35 kg
User length, approx.:	70-110 cm	90-115 cm	110-130 cm
Seat width:	25-28 cm	28-31 cm	31-34 cm
Seat depth:	24-27.5 cm	27.5- 31,5 cm	31,5- 35.5 cm
Back height, to arm axilla:	20-30 cm	23- 35 cm	22- 42 cm
Axilla to back top:	12 cm	12 cm	14 cm
Hip angle:	-20° +20°	-0° +30°	-0° +30°



### Leg length difference


size 1 - 4,5cm  
size 2 - 4,5 cm  
size 3 - 4,5 cm



### Abduction

size 1 - 40°  
size 2 - 35°  
size 3 - 25°


### HEAD SUPPORT



Dimensions in cm

	A	B
size 1	22	15
size 1 long	28	16
size 2	24	19
size 2 long	28	21


### BACK SUPPORT



Dimensions in cm

	A1	A2	B	C
size 1	22	19	31	10
size 2	27	24	35	11
size 3	29	26	40	13

### FOOT SUPPORT



Dimensions in cm

	A	B	C	D
size 1	31	19	15,5	19,5
size 2	34	22	17,5	22
size 3	39,5	25,5	19,5	26