




The Starfish is a development of the Zitzi Clozitt series. The seat has the same basic functions, but additional adjustment options. The seat can be adjusted to accommodate independent hip abduction or adduction as well as specific leg length discrepancies. There is also the option of catheterisation through a free space between the leg sections.


Catheterisation



Independent hip abduction or adduction as well as leg length discrepancy.

Adjustable tilt and hip angle

Active sitting and resting position!



Adjustable back

There is also the option of an adjustable back with thoracic pieces that are independently adjustable in both height and width. This option is particularly suitable for clients with a scoliosis.

