

# Zitzi Clozitt Travel



## Portable toilet and hygiene seat

Zitzi Clozitt Travel is a light and flexible toilet & hygiene seat, suitable for clients that require portability with specific positioning options.

## When you need a hygiene seat that

- Is simple and easy to transport
- Is easy to assemble and disassemble
- Is easy to mount on the toilet or the bathtub.

*Inspired by children and young people  
Developed by Anatomic SITT*

A Travel with "bathtub holders"



## Features of the Zitzi Travel

### Head support:

Padded head support which is height-adjustable.

### Back section:

The back section is also height-adjustable and is available in two sizes.

### Positioning support:

Keeps the child safely in position in the seat

### Seat section:

The seat section provides maximum pressure relief and numerous adjustment options

### Mounting:

Easily and safely mounted on the toilet

### Padding:

The padding is removable for hygiene reasons

### Portable:

Easily transportable and assembled



## Technical data & Article information

User length, approx.:	90-130 cm
Seat width:	28-34 cm
Seat depth:	27,5- 34,5 cm
Body width: Back size 2	27 cm
Body width: Back size 3	29 cm
Back height, to axil:	23- 36 cm
Abduction:	25°



We are certified according to ISO14001, Eco Management System

Box 6137, SE - 600 06 NORRKÖPING, Phone +46 (0) 11-16 18 00, fax +46 (0)11-16 20 05

[www.anatomicstt.com](http://www.anatomicstt.com)

**Anatomic**  
**-SITT-**