


### *"Flex-back"*

The back of the carseat rests against the seat of the car, and can be angled to provide an active sitting position as well as a resting position.



### Height-adjustable back

The height of the back can be easily adjusted so that the seat will grow along with the child.



### Adjustable leg sections

The leg pieces can be adjusted to accommodate independent hip abduction or adduction, as well as accommodate specific leg length discrepancies. An adjustment can also be made to the leg pieces to accommodate a pelvic obliquity. With these adjustments we can create an optimum seating position and this allows the chair to grow with the child.



Abduction



Leg length difference