


TECHNISCHE DATEN

Zitzi Carseat

User length, approx.:	125- 160 cm
Seat width:	33,5 cm
Seat depth, size	39 cm
Seat depth:	39- 45 cm
Body width:	29 cm
Back height, to axille:	39-45 cm
Hip angle:	-10° +40°
User weight:	15 - 36 kg
Weight:	8 kg


Zitzi Carseat Mini

	size 2	size 3
User length, approx.:	90-115 cm	110-130 cm
Seat width:	33,5 cm	33,5 cm
Seat depth, size	31 cm	35 cm
Seat depth:	29-33,5 cm	33,5-38cm
Body width:	27 cm	29 cm
Back height, to axille:	23- 35 cm	22- 42 cm
Hip angle:	-10° +40°	-10° +40°
User weight	9- 36 kg	15- 36 kg




HEAD SUPPORT

Carseat mini




Dimensions in cm	A	B
Carseat mini	27	20
Carseat	31	20

BACK SUPPORT



Dimensions in cm	A	B	C
size 1	22	31	10
size 2	27	35	11
size 3	29	40	13

FOOT SUPPORT



Dimensions in cm	A	B
size 1	34	22