


TECHNICAL DATA


This is how we have measured:

- 1 = Diagonally
- 2 = Width
- 3 = Height


Sternum support	1	2	3
X- Small:	23	7	10 cm
Small:	33	7	14 cm
Medium:	39	7	17 cm
Large:	51	7	20 cm




Sternum vest	1	2	3
X- Small:	23	7	10 cm
Small:	33	7	14 cm
Medium:	39	7	17 cm
Large:	51	7	20 cm



Cross strap	1	2
Small:	160	50 mm
Large:	250	70 mm



H support	1	2
X-small:	19,5	15-27 cm
Small:	26	17-31 cm
Medium:	34	18-38 cm
Large:	41	20-40 cm



	Measures around instep:
Instep holder, Small	140 - 170 mm
Instep holder, Medium	170 - 200 mm
Instep holder, Large	200 - 250 mm

	Measures around Bracelet:
Bracelet, Adult / Large	200 - 50 mm
Bracelet, Child / Small	160 - 50 mm

	Measurements, length of hand:
Hand tourniquet, Small	0 - 100 mm
Hand tourniquet, Medium	100 - 150 mm
Hand tourniquet, Large	160 - 220 mm

SittBelt, Adult / Large - S742	160 - 70 mm
SittBelt, Adult / Large - Car lock	160 - 70 mm
SittBelt, Child / Small . Fix lock	120 - 50 mm



We are certified according to ISO14001, Eco Management System

Box 6137, SE - 600 06 NORRKÖPING, Phone +46 (0) 11-16 18 00, fax +46 (0)11-16 20 05

www.anatomicstt.com

